

## F 93 X-ray Tank Fixer

### Description :

- X - ray fixer for tank processing of both medical and technical X - ray films;
- powdery;
- one - part.

### Processing Conditions :




fixing temperature :	18 °C	fixing time :	5 - 6 min.
	20 °C		4 - 5 min.
	22 °C		3 - 4 min.
	24 °C		2 - 3 min.

### Replenishment :

- For processing of 1 m<sup>2</sup> X - ray film **about 600 ml** tank solution are necessary at standard conditions.

### Mix Instructions :

preparation X - ray fixer :

Water ( 30 - 40 °C )	F 93 - Powder	fill up with Water to
 about ¾ of tank - , replenisher volume resp.		
7,5 l	Dissolve the powder under permanent stirring and then fill up to the final volume and mix good .	10,0 l
17,0 l		22,5 l
11,0 l		15,0 l

**F 93 X-ray Tank Fixer**

**Assortment :**

Order Number	Product Name	Mix Size for	Capacity ( m <sup>2</sup> )
10193	F 93	10 x 10,0l	165
10194		5 x 22,5l	185
10195		7 x 15,0l	175

**Storage and Handling :**

- Until the preparation F 93 has to be kept in the closed original package and stored at temperatures of 10 to 25 ° C
- All users of photographic processing solutions should exercise the greatest care to avoid the chemicals contacting the skin or eyes . Always wear solution resistant gloves .
- The disposal of used photographic processing solutions has to be occurred in accordance with the legal regulations .