

## T 93 X-ray Tank Developer

### Description :

- developer for (manual) tank processing of both medical and technical X - ray films;
- powdery;
- two - part.

### Processing Conditions :





developing temperature :	18 °C	developing time :	5 - 6 min.
	20 °C		4 - 5 min.
	22 °C		3 - 4 min.
	24 °C		2 - 3 min.

### Replenishment :

- For processing of 1 m<sup>2</sup> X - ray film 300 to 400 ml developer solution are necessary at standard conditions.

### Mix Instructions :

preparation of developer :

Water (30 - 40 °C)	Part 1	Part 2	fill up with Water to final volume of ...
 about ¾ of tank - , replenisher volume resp.			
7,5 l	First dissolve part 1 completely in water , afterwards add part 2 under permanent stirring, then fill up to the final volume and mix good.		10,0 l
17,0 l			22,5 l
11,0			15,0 l

**T 93 X-ray Tank Developer**

**Assortment :**

Order Number	Product Name	Mix Size For	Capacity ( m <sup>2</sup> )
10133	T 93	10 x 10,0l	250 - 330
10134		5 x 22,5l	280 - 375
10135		7 x 15,0l	260 - 350

**Storage and Handling :**

- Until the preparation T 93 has to be kept in the closed original package and stored at temperatures of 10 to 25 ° C ..
- All users of photographic processing solutions should exercise the greatest care to avoid the chemicals contacting the skin or eyes . Always wear solution resistant gloves .
- The disposal of used photographic processing solutions has to be occurred in accordance with the legal regulations .