

**UNIROLL-D X-ray Developer for Roller Machines**

**Description :**






- X - ray developer for all customary developing machines with roller transport and for spray-developing machines ;
- good sensitometric results with low fog values in a wide range of temperatures;
- Because of the special oxidation stability the developer is especially suitable for machines with low and / or irregular utilization.
- liquid concentrate;
- three - part;

**Replenishment :**

- For processing of 1 m<sup>2</sup> X-ray film 400 to 600 ml replenisher are necessary at standard conditions.
- Generally the replenishment occurs automatically by addition of the solution in dependence of the processing film.

**Mix Instructions :**

preparation X - ray developer replenisher :

Water ( 25 - 30 °C )	Part 1	Part 2	Part 2	fill up with Water to
				
14 l	5,0 l	0,5 l	0,5 l	20 l
7 l	2,5 l	0,25 l	0,25 l	10 l

First dissolve part 1 completely in water, afterwards add part 2 under permanent stirring and then part 3 in the same way.

preparation X - ray developer working solution :

20 ml Uni - Starter per litre contents of tank has to be added the prepared replenisher .

## UNIROLL-D X-ray Developer for Roller Machines

### Assortment :

Order Number	Product Name	Mix Size for	Concentrate	Capacity m <sup>2</sup>
10510	UNIROLL-D	2 x 20 l	Teil 1 2 x 5,0 l Teil 2 2 x 0,5 l Teil 3 2 x 0,5 l	65 - 100
10511	UNIROLL-D	3 x 10 l	Teil 1 3 x 2,5 l Teil 2 3 x 0,25l Teil 3 3 x 0,25l	50 - 75

### Storage and Handling :

- Until the preparation the concentrate has to be kept in the closed original package and stored at temperatures of 10 to 25 °C .
- **UNIROLL-D** contains components which can cause irritations in case of contact with the skin.
- In case of contact with the eyes it has to be rinsed thoroughly with plenty of running cold water by opened eyelid and a medical doctor has to be consulted.