

UNIROLL-F X-ray Fixer for Roller Machines

Description :





- X - ray fixer for all customary developing machines with roller transport and for spray - developing machines with a fixage time by 15 s or more ;
- Because of ist excellent hardening properties and the fast working fixer components the product is suitable for use in fast going machines .
- application in conjunction with Developers, like XROLL-D and UNIROLL-D
- liquid concentrate ;
- two - part ;

Replenishment :

- For fixing of 1 m² X - ray film **600 to 750 ml** replenisher are necessary at standard conditions .
- Generally the replenishment occurs automatically by addition of the solution in dependence of the processing film .

Mix Instructions :

preparation X - ray fixer replenisher :

Water (25 - 30 °C)	Part 1	Part 2	fill up with Water to
			
14,0 l	5,0 l	1,0 l	20,0 l
7,0 l	2,5 l	0,5 l	10,0 l

- First dissolve part 1 completely in water, afterwards add part 2 under permanent stirring ;
- When using an adapter **UNIROLL-F** is also suitable for Agfa- and Kodak - Mixer .

UNIROLL-F X-ray Fixer for Roller Machines

Assortment :

Order Number	Product Name	Mix Size for	Concentrate	Capacity m ²
10570	UNIROLL-F	2 x 20 l	Teil 1 2 x 5,0 l Teil 2 2 x 1,0 l	50 - 65
10571		3 x 10 l	Teil 1 3 x 2,5 l Teil 2 3 x 0,5 l	40 - 50

Storage and Handling :

- Until the preparation the concentrates has to be kept in the closed original package and stored at temperatures of 10 to 25 °C .
- **UNIROLL-F** contains components which can cause irritations in case of contact with the skin.
- In case of contact with the eyes it has to be rinsed thoroughly with plenty of running cold water by opened eyelid and a medical doctor has to be consulted.