

Fixing Bath and Replenisher C-41 RA FX-R

Description :




- fixing bath for the process C-41 RA / RANP;
- low fixing time;
- one-part;

Technical Data :

Processing Time:	1 min 30 s
Temperature:	35 -41 °C
Replenishment Rate:	33 ml / film 135-24

Mix Instructions:

preparation fixing bath - replenisher:

Water (20-30 °C)	C-41 RA FX-R Conc.	fill up with Water to
		
1,5 l	2,0 l	4 l
2,0 l	2,5 l	5 l
4,5 l	5,0 l	10 l

preparation fixing bath - working solution:

1000 ml replenisher = 1000 ml working solution

Fixing Bath and Replenisher C-41 RA FX-R

Assortment :

Order Number	Product Name	Mix Size for	Concentrate	Capacity Films 135-24
13350	C-41 RA FX-R	2 x 10 l	2 x 5,0 l	606
13351		2 x 5 l	2 x 2,5 l	303
13353		5 x 4 l	5 x 2,0 l	606