

Stabilizer and Replenisher RA-4 STAB-WL

Description :

- stabilizer bath for the washless RA-4-process ; increase of the colour stability ;
- formaldehyde free ;
- concentrate in a double-necked bottle with dosing chamber, mixing doses to make 5 L or 10 L of replenisher,
- whole bottle content is to make 20 x 5 L / 10 x 10 L or 100 L replenisher

Technical Data :

Processing Time :	1 min 30 s	Temperature :	29 - 38 °C
Replenishment Rate :		250 ml / m ²	at 4 tanks
		360 ml / m ²	at 3 tanks

Mix instruction / Dosage:

	5 l	10 l	100 l
1. Start with water (20-40 °C)			
2. Open cap of dosing chamber			
Pump up 25 mls of concentrate into dosing chamber by pressing the bottle (see figure...)	... 1 x	... 2 x	
4. Empty dosing Chamber	... 5 L	... 10 L	Take whole content
5. Mix			

preparation stabilizer - working solution :

1000 ml r e p l e n i s h e r = 1000 ml w o r k i n g s o l u t i o n

Stabilizer and Replenisher RA-4 STAB-WL

Assortment :

REF # Number	Product Name	Mix Size to make ... per pack	Contents of Concentrate per pack	Capacity (m ²)
16033	RA-4 STAB - WL	2 x (20 x 5 l) or 2 x (10 x 10 L) or 2 x 100 l	2 x 0,5 l	554 - 800